# API 20 E

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- 35 ~37 °C incubator
- Refrigerator
- Bunsen burner
- Marker pen

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Suspension Medium 5 ml	Demineralized water	
TDA reagent 5 ml	Ferric chloride H <sub>2</sub> O	3.4 g qsp 100 ml
JAMES reagent 5 ml	Compound J 2183 (confidential) HCl 1N	0.5 g qsp 100 ml
VP 1 reagent 5 ml	Potassium hydroxide H <sub>2</sub> O	40 g 100 ml
VP 2 reagent 5 ml	α-naphthol Ethanol	6 g 100 ml
NIT 1 reagent 5 ml	Sulfanilic acid Acetic acid H <sub>2</sub> O	0.4 g 30 g 70 ml
NIT 2 reagent 5 ml	N,N-dimethyl-1-naphthylamine Acetic acid H <sub>2</sub> O	0.6 g 30 g 70ml
Zn reagent	Zinc dust	10g

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TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS / ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
ONPG	2-nitrophenyl-βD-galactopyranoside	0.223	β-galactosidase (Ortho NitroPhenyl-βD-Galactopyranosidase)	colorless	yellow (1)
<u>ADH</u>	L - arginine	1.9	Arginine DiHydrolase	yellow	red / orange (2)
<u>LDC</u>	L - lysine	1.9	Lysine DeCarboxylase	yellow	red / orange (2)
<u>ODC</u>	L - ornithine	1.9	Ornithine DeCarboxylase	yellow	red / orange (2)
<u>CIT</u>	trisodium citrate	0.756	CITrate utilization	pale green / yellow	blue-green / blue (3)
<u>H<sub>2</sub>S</u>	sodium thiosulfate	0.075	H <sub>2</sub> S production	colorless / greyish	black deposit / thin line
<u>URE</u>	urea	0.76	UREase	yellow	red / orange (2)
TDA	L- tryptophane	0.38	Tryptophane DeAminase	TDA / immediate yellow   reddish brown	
IND	L- tryptophane	0.19	INDole production	JAMES / immediate colorless   pink pale green / yellow	
<u>VP</u>	sodium pyruvate	1.9	acetoin production (Voges Proskauer)	VP 1 + VP 2 / 10 min colorless   pink / red (5)	
<u>GEL</u>	Gelatin (bovine origin)	0.6	GELatinase	no diffusion	diffusion of black pigment
GLU	D-glucose	1.9	fermentation / oxidation (GLUcose) (4)	blue / blue-green	yellow / greyish yellow
MAN	D-mannitol	1.9	fermentation / oxidation (MANnitol) (4)	blue / blue-green	yellow
INO	inositol	1.9	fermentation / oxidation (INOsitol) (4)	blue / blue-green	yellow
SOR	D-sorbitol	1.9	fermentation / oxidation (SORbitol) (4)	blue / blue-green	yellow
RHA	L-rhamnose	1.9	fermentation / oxidation (RHAmnose) (4)	blue / blue-green	yellow
SAC	D-sucrose	1.9	fermentation / oxidation (SACcharose) (4)	blue / blue-green	yellow
MEL	D-melibiose	1.9	fermentation / oxidation (MELibiose) (4)	blue / blue-green	yellow
AMY	amygdalin	0.57	fermentation / oxidation (AMYgdalin) (4)	blue / blue-green	yellow
ARA	L-arabinose	1.9	fermentation / oxidation (ARAbinose) (4)	blue / blue-green	yellow
OX	(see oxidase test package insert)		cytochrome - OXidase	(see oxidase test package insert)	

## SUPPLEMENTARY TESTS

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
Nitrate reduction GLU tube	potassium nitrate	0.076	NO <sub>2</sub> production reduction to N <sub>2</sub> gas	NIT 1 + NIT 2/2-5 min yellow   red	
				orange-red	Zn / 5 min yellow
MOB	API M Medium or microscope	-	motility	non-motile	motile
McC	MacConkey medium	-	growth	absence	presence
OF-F OF-O	glucose (API OF Medium)	-	fementation : under mineral oil oxidation : exposed to the air	green green	yellow yellow

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